

RealLine Ureaplasma urealyticum / Ureaplasma parvum Str-Format

Instructions for Use

RealLine Ureaplasma urealyticum / U. parvum Str-Format

A QUALITATIVE ASSAY KIT FOR THE DIFFERENTIAL DETECTION OF *UREAPLASMA UREALYTICUM* DNA AND *UREAPLASMA PARVUM* DNA BY REAL TIME PCR METHOD










In vitro Diagnostics



RealLine Ureaplasma urealyticum / U. parvum (Str-Format)	VBD2294	96 Tests
valid from	September 2019	

RealLine Ureaplasma urealyticum / U. parvum Str-Format

Explanation of symbols used in labeling

	<i>In vitro</i> diagnostic medical device
	Batch code
	Catalogue number
	Contains sufficient for <n> tests
	Use-by-date
	Temperature limit
	Consult instructions for use
	Keep away from sunlight
	Manufacturer



BIORON Diagnostics GmbH

In den Rauhweiden 20
67354 Römerberg
Germany

Phone +49 6232 298 44 0
Fax: +49 6232 298 44 29
info@bioron.de

RealLine *Ureaplasma urealyticum* / *U. parvum* Str-Format

Table of content:

1. INTRODUCTION	4
2. KIT CONTENTS	5
3. PRINCIPLE OF THE METHOD	5
4. SPECIFICATIONS	6
5. WARNING AND PRECAUTIONS	7
6. ADDITIONAL MATERIALS AND DEVICES REQUIRED BUT NOT SUPPLIED	7
7. PREPARATION OF SPECIMENS	8
8. PROCEDURE	8
9. DATA ANALYSIS AND INTERPRETATION	10
10. STORAGE AND TRANSPORTATION	11
ANNEX I: Settings for RealLine Cyclor and DT96:	11

RealLine *Ureaplasma urealyticum* / *U. parvum* Str-Format

A QUALITATIVE ASSAY KIT FOR THE DIFFERENTIAL DETECTION OF *UREAPLASMA UREALYTICUM* DNA AND *UREAPLASMA PARVUM* DNA BY REAL TIME PCR METHOD

In vitro diagnosticum

1. INTRODUCTION

Clinical Information:

Ureaplasma urealyticum is part of the normal genital flora of men and women and is found in about 70 % of sexually active humans. *U.urealyticum* can be associated with a number of diseases in humans, including non-specific urethritis (NSU), infertility, chorioamnionitis, stillbirth, premature birth, and, in the perinatal period, pneumonia, bronchopulmonary dysplasia and meningitis.

Ureaplasma parvum colonized in women up to 57 % in the vagina and can result during a pregnancy - caused by inflammation - to premature rupture of membranes, premature labor, chorioamnionitis, and the child to sepsis and bronchopulmonary dysplasia and to meningitis.

RealLine *Ureaplasma urealyticum* / *Ureaplasma parvum* assay kit is intended for the differential detection of *Ureaplasma urealyticum* and *Ureaplasma parvum* DNA in clinical specimens (urogenital and cervical swabs, semen, prostate fluid, urine) using the method of real-time polymerase chain reaction (PCR) with fluorescence detection of amplified product.

RealLine *Ureaplasma urealyticum* / *U. parvum* (Str-format) assay kit is intended for detection of *Ureaplasma urealyticum* DNA and *Ureaplasma parvum* DNA isolated from specimens using extraction kits:

RealLine DNA-Express (REF VBC8899)

RealLine DNA-Extraction 2 (REF VBC8897)

RealLine DNA-Extraction 3 (REF VBC8889)

RealLine Extraction 100 (REF VBC8896)

The **Str-Format Kit** contains 96 tubes (0.2 ml) in strips with lyophilized Mastermix. 50 µl of extracted DNA have to be pipetted into the tube and the ready mastermix is diluted. The kit contains reagents required for 96 tests, including control samples and the positive control sample.

The kit is intended to use with block cyclers: iQ™5 iCycler, CFX™96 (Bio-Rad, USA), DT-96 (DNA-Technology, Russia) and RealLine Cycler (BIORON Diagnostics GmbH).

The use of:

- ! **Extraction Kits for nucleic acids from clinical specimen from other supplier**
- ! **other real-time PCR devices**
- ! **appropriate reaction volumes, other than 50 µl**

has to be validated in the lab by the user. The special notes regarding the internal control IC have to be strongly followed.

RealLine Ureaplasma urealyticum / U. parvum Str-Format

2. KIT CONTENTS

Universal Positive Control Sample (PC)	1 vial, 1 ml
Ready Master Mix (RMM) , lyophilized	96 test-tubes 12 strips x 8 tubes
The kit is additionally supplied with optical-quality PCR-film	

3. PRINCIPLE OF THE METHOD

The Real-time PCR is based on the detection of the fluorescence, produced by a reporter molecule, which increases as the reaction proceeds. Reporter molecule is dual-labeled DNA-probe, which specifically binds to the target region of pathogen DNA. Fluorescent signal increases due to the fluorescent dye and quencher separating by Taq DNA-polymerase exonuclease activity during amplification. PCR process consists of repeated cycles: temperature denaturation of DNA, primer annealing and complementary chain synthesis.

Threshold cycle value – Ct – is the cycle number at which the fluorescence generated within a reaction crosses the fluorescence threshold, a fluorescent signal rises significantly above the background fluorescence. Ct depends on initial quantity of pathogen DNA template.

The use of **Internal Control (IC)** prevents generation of false negative results associated with possible loss of DNA template during specimen preparation. IC indicates if PCR inhibitors occur in the reaction mixture. IC template should be added in each single sample (including control samples) prior to DNA extraction procedure. The amplification and detection of IC does not influence the sensitivity or specificity of the target DNA PCR.

Note: IC is a component of the NA extraction kits of RealLine series. Internal Control is added to the sample during NA isolation step and is used throughout the whole process of NA extraction, amplification, detection.

RealLine *Ureaplasma urealyticum* / *U. parvum* Str-Format

4. SPECIFICATIONS

4.1 Sensitivity:

Sensitivity of the detection of *Ureaplasma urealyticum* DNA is determined on five samples prepared from *Ureaplasma urealyticum* DNA SRS (Standard Reference Sample), containing 100 copies of *Ureaplasma urealyticum* DNA per sample. Sensitivity equals 100%.

Sensitivity of *Ureaplasma parvum* DNA detection is determined on five samples prepared from *Ureaplasma parvum* DNA SRS, containing 100 copies of *Ureaplasma parvum* DNA per sample. Sensitivity equals 100%.

4.2 Specificity:

Specificity of *Ureaplasma urealyticum* DNA detection equals 100 %. It is determined on five samples prepared from *Ureaplasma parvum* DNA SRS, containing 100 copies of *Ureaplasma parvum* DNA per sample.

Specificity of *Ureaplasma parvum* DNA detection equals 100 %. It is determined on five samples prepared from *Ureaplasma urealyticum* DNA SRS, containing 100 copies of *Ureaplasma urealyticum* DNA per sample.

4.3 Diagnostic sensitivity:

Diagnostic sensitivity of *Ureaplasma urealyticum* DNA detection: clinical trials conducted on 47 positive samples showed 100 % sensitivity (interval 93.8 % -100 %, with a confidence level of 90 %);

Diagnostic sensitivity of *Ureaplasma parvum* DNA detection: clinical trials conducted on 75 positive samples showed 100% sensitivity (interval 96.1 % -100 %, with a confidence level of 90 %);

4.4 Diagnostic specificity

Diagnostic specificity of *Ureaplasma urealyticum* DNA detection: clinical trials conducted on 130 negative samples showed 100 % specificity (97.7 % -100 % interval, with a 90 % confidence level).

Diagnostic specificity of *Ureaplasma parvum* DNA detection: clinical trials conducted on 113 negative samples showed 100% specificity (interval 97.4 % -100 %, with a confidence level of 90 %).

Analysis by the CE-marked reference kit showed full match of results.

RealLine Ureaplasma urealyticum / U. parvum Str-Format

5. WARNING AND PRECAUTIONS

- ☞ For In vitro use only.
- ☞ The kits must be used by skilled personnel only.
- ☞ When handling the kit, follow the national safety requirements for working with pathogens.
- ☞ To prevent contamination, the steps of DNA isolation and PCR test run must be spatially separated.
- ☞ Avoid microbial and ribonuclease contamination of reagents when removing aliquots from reagent vials.
- ☞ Wear protective disposable gloves, laboratory coats and eye protection when handling specimens and kit reagents.
- ☞ Every workplace must be provided with its own set of variable-volume pipettes, necessary auxiliary materials and equipment. It is prohibited to relocate them to other workplaces.
- ☞ The use of sterile disposable pipette tips is recommended.
- ☞ Never use the same tips for different samples.
- ☞ Do not pool reagents from different lots or from different vials of the same lot.
- ☞ Dispose unused reagents and waste in accordance with country, federal, state and local regulations.
- ☞ Do not use the kit after the expiration date at the side label of the kit.

6. ADDITIONAL MATERIALS AND DEVICES REQUIRED BUT NOT SUPPLIED

- real time PCR system, like described in p.1
- DNA-Extraction Kit: **RealLine DNA-Express**, **RealLine DNA-Extraction 2**, **RealLine DNA-Extraction 3** or **RealLine Extraction 100**
- Internal Control reagent (VBC8881) and Negative Control Sample, if the kit is used with the extraction kits of other supplier.
- laminar safety box;
- refrigerator;
- half-automatic variable-volume single-channel pipettes;
- disposable medical non-sterile powder-free gloves;
- disposable pipette tips with aerosol barrier;
- biohazard waste container
- razor or scalpel.

RealLine Ureaplasma urealyticum / U. parvum Str-Format

7. PREPARATION OF SPECIMENS

Each group of samples undergoing the procedure of DNA isolation must include a **Positive Control sample (PC)** from this kit and a **Negative Control sample (NC)** which is a component of the DNA extraction kit.

We strongly recommend the implementation of the Internal Control IC, the Negative Control NC and Positive Control PC samples to the extraction procedure.

When using a kit of another supplier for the extraction of nucleic acids as recommended in p.1. add **20 µl** of **IC (VBC8881)** to each tube.

- For the NC use **100 µl** of Negative Control Sample.
- For the PC use **70 µl** of Negative Control Sample and **30 µl** of Positive Control to the tube marked PC.

The assay is performed on extracted DNA specimens obtained from the clinical material using one of the DNA extraction kits listed in p.1.

8. PROCEDURE

8.1. Preparation of the kit components:

Prior to the test take the kit out of the refrigerator and keep the **Ready Master Mix (RMM)** closed in the package at (18 – 25) °C for at least 30 minutes. Then open the package and cut the necessary number of tubes in strips with RMM (*for prepared samples and controls: 1 PC and 1 NC*) with the razor or scalpel. Cut the tubes together with the covering film.

Put the remaining strips immediately back into the foil pouch, squeeze the air out and tightly close with the clip.

After initial opening store RMM at (2 – 8) °C for no more than 3 months.

After initial opening, store PC at (2 - 8) °C for no more than 1 month or in 50 µl aliquots at (-18 ...-24) °C for no more than 3 months.

8.2. Label the tubes with RMM for each specimen and control.

Attention! Labels should be placed on the lateral side of the tubes.

8.3. Add **50 µl** of corresponding extracted DNA solution to each tube using a separate pipette tip with filter. Tightly close the tubes with caps or seal with the PCR optical quality film.

8.4. Place the tubes into the Real Time PCR system.

RealLine *Ureaplasma urealyticum* / *U. parvum* Str-Format

8.5. Program Real Time PCR system as follows:

Step 1:	50°C	2min	50 cycles
Step 2:	95°C	2min	
Step 3:	94°C	10 sec	
	60°C*	20 sec	

* Measure the fluorescence at 60°C

8.6. Select the amplification detection channels:

- Collect Real Time PCR data through the **FAM** channel for detection of amplification of **IC DNA**.
- Collect Real Time PCR data through the **HEX** channel for detection of amplification of ***Ureaplasma urealyticum* DNA**.
- Collect Real Time PCR data through the **ROX** channel for detection of amplification of ***Ureaplasma parvum* DNA**.

8.7. Program the positions of test tubes with samples, positive and negative controls according to the instruction manual for the Real Time PCR system in use.

8.8. Run the program.

RealLine *Ureaplasma urealyticum* / *U. parvum* Str-Format

9. DATA ANALYSIS AND INTERPRETATION

9.1. For **PC** the program should detect:

- increase of the IC DNA amplification signal (channel **FAM**) and determine the IC Ct;
- increase of the *Ureaplasma urealyticum* DNA amplification signal (channel **HEX**) and determine the PC Ct value.
- increase of the *Ureaplasma parvum* DNA amplification signal (channel **ROX**) and determine the PC Ct value;

9.2. For **NC** the program should detect the increase of the amplification signal of IC DNA (channel **FAM**) and determine the threshold cycle, IC Ct. No **ROX** or **HEX** fluorescent increase should appear.

When **Ct** value for NC through **HEX and ROX** channel is \leq **Ct 40**, this indicates the presence of contamination (see paragraph 9.8).

9.3. For each sample the program should detect the increase of the amplification signal of IC DNA (channel **FAM**) and determine IC **Ct**.

9.4. Calculate (IC **Ct**)_{av} as an average IC **Ct** of all analyzed samples (including PC and NC). IC **Ct** values that differ by more than 2 from the (IC **Ct**)_{av} should be ignored. Recalculate the (IC **Ct**)_{av} for the remaining values after the screening.

9.5. The specimen is considered **negative** (not containing *Ureaplasma urealyticum* / *Ureaplasma parvum* DNA), if **Ct** value via **HEX** and **ROX** channel for this sample is **above Ct 40** or not determined.

If IC **Ct** value for such sample differs from the (IC **Ct**)_{av} value by more than 2, the result is regarded as equivocal. A repeated analysis of the sample, starting with the DNA isolation step is necessary.

9.6. If **Ct** value via **HEX** channel for this sample is **less than or equals to Ct 40**, the sample is considered **positive** and contains *Ureaplasma urealyticum* DNA,

If **Ct** value via **ROX** channel for this sample is **less than or equals to Ct 32**, the sample is containing *Ureaplasma parvum* DNA and considered **positive**.

If **Ct** value via **ROX** channel for this sample is **> Ct 32 and \leq Ct 40**, the sample is containing *Ureaplasma parvum* DNA and considered as **positive**, but possible clinically insignificant.

Note: *Ureaplasma parvum* is a commensal bacterium in the uterus as part of the microbiome in healthy women. In low concentrations these bacteria do not necessarily need a treatment.

The **RealLine *Ureaplasma urealyticum* / *Ureaplasma parvum*** kit detects even lowest amounts of *Ureaplasma parvum*. For results with **Ct 32 to 40** it is advisable to re-check the samples regarding the subsequent treatment decisions, by taking the clinical data, course and appearance of the patient, into account.

RealLine Ureaplasma urealyticum / U. parvum Str-Format

9.7. If **Ct** value for NC through **HEX and ROX** channel is \leq **Ct 40**, this indicates the presence of contamination. In case of contamination all positive results of this individual PCR run are considered equivocal. Actions are required to identify and eliminate the source of contamination, and repeat the analysis of all samples of this run that were identified as positive. Samples that showed negative results in this run should be considered as negative.

10. STORAGE AND TRANSPORTATION

- Transport and store the assay kit at (2 – 8) °C in the manufacturer's packing.
- Transportation at up to 25 °C for up to 10 days is allowed.
- Do not freeze the kit!
- Do not pool reagents from different lots or from different vials of the same lot.
- Strictly follow the Instruction manual for reliable results.
- Do not use kits with damaged inner packages and get in contact with BIORON Diagnostics GmbH.

- **Storage and shelf life of solutions and components of the kit after initial opening:**
 Positive Control sample: 1 month at (2 – 8) °C, in 50 µl aliquots 3 month at (-18 ...-24)°C .
 Ready Master Mix (RMM): 3 months at (2 – 8) °C.

Technical support: techsupport@bioron.de

ANNEX I: Settings for RealLine Cyclers and DT96:

for these cyclers the measurement exposure must be adjusted. Choose the **Operation with the device** mode in the **Settings** menu, select the item **Measurement exposition:**

- **FAM to 250**
- **HEX and ROX to 1000**

Confirm that the current exposure value is saved by pressing **YES**

Attention! The specified exposure values are applicable only for RealLine kits and, if necessary, must be changed for other purposes.

**RealLine Ureaplasma urealyticum / U. parvum
Str-Format**

