

RealLine Pathogen Diagnostic Kits

for Real Time PCR

SEXUALLY TRANSMITTED INFECTIONS
HUMAN PAPILLOMA VIRUS
HUMAN HERPES VIRUS
... AND OTHERS



 **BIORON**
diagnostics



Your Partner for Real Time PCR Detection

BIORON Diagnostics GmbH is manufacturer and worldwide acting trader of diagnostic products.

RealLine kits are based on

- the most advanced variants of PCR
- fast and sensitive technique

BIORON products fulfil highest quality and health standards.

Our customers value us for our individual concern and prompt support.


Our QM-System has been certified by



RealLine Pathogen Kits

- Lyophilized Mastermixes
- Transport at room temperature (<10 days)
- Ready for use
- Simple procedure
- Same protocols for PCR and RT-PCR
- Controls included
- Less handling steps
- Less contamination risk
- CE-IVD

Fast and easy



Open System

- Usable with various real-time cyclers

Sexually Transmitted Infections

Single Kits	Str-format	Fla-format
RealLine Chlamydia trachomatis	VBD1998	VBD1996
RealLine Neisseria gonorrhoeae	VBD4498	VBD4496
RealLine Mycoplasma hominis	VBD4398	VBD4393
RealLine Mycoplasma genitalium	VBD4396	VBD4394
RealLine Ureaplasma urealyticum	VBD2298	VBD2296
RealLine Trichomonas vaginalis	VBD2098	VBD2096
RealLine Gardnerella vaginalis	VBD4598	VBD4596
RealLine Candida albicans	VBD4698	VBD4696
RealLine Treponema pallidum	VBD1898	VBD1896

Multiplex Kits	Str-format	Fla-format
RealLine Chlamydia trachomatis / Neisseria gonorrhoeae	VBD0457	VBD0458
RealLine Chlamydia trachomatis / Mycoplasma genitalium	VBD0490	VBD0487
RealLine Chlamydia trachomatis / Ureaplasma urealyticum	VBD0492	VBD0491
RealLine Mycoplasma genitalium / Mycoplasma hominis	VBD0494	VBD0493
RealLine Ureaplasma urealyticum / Ureaplasma parvum	VBD2294	VBD2295
RealLine Trichomonas vaginalis / Gardnerella vaginalis	VBD0477	VBD0478
RealLine Trichomonas vaginalis / Neisseria gonorrhoeae	VBD0496	VBD0495
RealLine Candida albicans / Gardnerella vaginalis	VBD0444	VBD0445



Human Papilloma Virus

High Carcinogenic Risk Types	Str-format	Fla-format
RealLine HPV 16 / 18	VBD8473	VBD8474
RealLine HPV 31 / 33	VBD8471	VBD8480
RealLine HPV 35 / 45	VBD8469	VBD8477
RealLine HPV HCR Genotype	VBD8479	VBD8482
RealLine HPV HCR Genotype quantitative	VBD8478	–

The genotyping kits are detecting the Human Papilloma Types: 16, 18, 31, 33, 35, 39, 45, 51, 52, 56, 58, 59

Low Carcinogenic Risk Types	Str-format	Fla-format
RealLine HPV 6 / 11	VBD8475	VBD8481

Human Herpes Virus

Product	Str-format	Fla-format
RealLine HSV1/ HSV2	VBD2195	VBD2197
RealLine CMV quantitative	VBD1598	VBD1596
RealLine EBV quantitative	VBD2198	VBD2196
RealLine HHV-6	VBD2150	VBD2151
RealLine VZV	VBD2185	VBD2187

Tick Borne Infections

Product	Str-format	Fla-format
RealLine TBEV (FSME)	VBD1198	VBD1199
RealLine Borrelia burgdorferi s.l.	VBD1498	VBD1499
RealLine Anaplasma / Ehrlichia	VBD5398	VBD5399

... and others

Product	Uni-format	Str-format	Fla-format
RealLine Toxoplasma gondii	–	VBD1798	VBD1796
RealLine Rubella	–	VBD2598	VBD2596
RealLine Helicobacter pylori	–	VBD3798	VBD3796
RealLine Mycobacterium tuberculosis	VBD2398*	VBD2399	–

*Including extraction

RealLine Pathogen Extraction Kits

Product	Clinical Material	MagicMag Rack	Nucleic Acid	Ref .No
RealLine DNA Express	epithelial cell specimen, urine, sperm	–	DNA	VBC8899
RealLine DNA Extraction 2	epithelial cell specimen, urine, liquor, biopsy material, blood serum, plasma	+	DNA	VBC8897
RealLine DNA Extraction 3	epithelial cell specimen, urine, sperm, blood serum, plasma	+	DNA	VBC8889
RealLine Extraction 100	epithelial cell specimen, liquor, biopsy material, blood serum, plasma, whole blood, ticks, urine, sperm, stool	+	DNA and RNA	VBC8896

Internal Control included in the Extraction Kits, but also available separately

Accessories

Product	Ref .No
RealLine Internal Control	VBC8881
MagicMag Rack	BI87000

Unique - available from BIORON only:

The tubes remain inside the rack for the whole purification process. This is possible due to the built-in tilt mechanism.



To facilitate your order,

please provide us with the following information

- Institution or company name
- Customer name, email address, telephone number
- Billing address
- Shipping address
- Valid purchase order number
- Product number and required quantity

Terms and conditions of BIORON Diagnostics GmbH

Unless agreed otherwise the single transactions between client and BIORON are to be carried out according to

- **International:** The General Conditions of Export
<http://www.bioron.de/index.php/en/gtc/gce-international/>
- **Germany:** Allgemeine Geschäfts- und Lieferbedingungen
<http://www.bioron.de/index.php/en/gtc/agb-national/>

Find our products in your country

BIORON Diagnostics is represented in various countries worldwide.

Please view the information regarding our distributors in your local area via


<http://www.bioron.de/index.php/en/distribution/>

If you do not find a distributor listed in your country please contact us directly in Germany.

Imprint

BIORON Diagnostics GmbH
Rheinhorststr. 18
67071 Ludwigshafen, Germany
Phone +49 621 545 900 70
Fax +49 621 545 900 68
info@bioron.de
www.bioron.de

Legal form:	Limited Liability Company
Managing Director:	Dr. Ferdinand Holzinger
General Manager / Prokurist:	Barbara Hassib
Commercial Register:	AG Ludwigshafen HRB 02298
Registered Office:	Ludwigshafen
Value Added Tax Identification Number: (VAT ID)	DE 274 094 961
German Tax Number:	27/656/01009



BIORON Diagnostics GmbH

Rheinhorststr. 18
67071 Ludwigshafen
Germany

Phone +49 621 545 900 70

Fax +49 621 545 900 68

info@bioron.de

www.bioron.de