

Human Papilloma Virus

HPV



BIORON Diagnostics offers a wide range of detection kits for different types of the Human papilloma Viruses i.e.

RealLine HPV HCR Screen Kit

which detects 14 different high carcinogenic risk types of the Human Papilloma Viruses in one test.

Types: 16, 18, 31, 33, 35, 39, 45, 51, 52, 56, 58, 59, 66, 68

Advantages of the RealLine HPV kits:

- Same Extraction and PCR Protocol
- Same internal Control System
- All CE-IVD marked



 **BIORON**
diagnostics



Human Papilloma Virus

High Carcinogenic Risk Types	Str-format	Fla-format
RealLine HPV 16 / 18	VBD8473	VBD8474
RealLine HPV 31 / 33	VBD8471	VBD8480
RealLine HPV 35 / 45	VBD8469	VBD8477
RealLine HPV HCR Screen	VBD8444	VBD8445
RealLine HPV HCR Genotype	VBD8479	VBD8482
RealLine HPV HCR Genotype quantitative	VBD8478	–

The genotyping kits are detecting the Human Papilloma Types: 16, 18, 31, 33, 35, 39, 45, 51, 52, 56, 58, 59

Low Carcinogenic Risk Types	Str-format	Fla-format
RealLine HPV 6 / 11	VBD8475	VBD8481

We recommend to use the following extraction:

Pathogen Extraction Kits	Clinical Material	Packsize	Nucleic Acid	Ref .No
RealLine DNA Express	epithelial cell specimen, urine, sperm	100 rcs	DNA	VBC8899
RealLine DNA Extraction 3	epithelial cell specimen, urine, sperm, blood serum, plasma	96 rcs	DNA	VBC8889
RealLine Extraction 100	epithelial cell specimen, liquor, biopsy material, blood serum, plasma, whole blood, ticks, urine, sperm, stool	48 rcs	DNA and RNA	VBC8896

Internal Control included in the Extraction Kits, but also available separately