

RealLine Extraction Kits – Overview

	Number of Reactions	REF	Clinical Material	Extraction of	Detection of	Magnet Rack
RealLine DNA Express CE-IVD	100	VBC8899	epithelial cell specimen, urine, sperm, saliva	DNA	STI, HHV, HPV, Parvovirus B19	no
RealLine DNA Extraction 2 CE-IVD	96	VBC8897	epithelial cell specimen urine, liquor, biopsy material, blood serum/plasma	DNA	STI, HHV, HPV, Parvovirus B19 Toxoplasma gondii Helicobacter pylori	yes
RealLine DNA Extraction 3 CE-IVD	96	VBC8889	epithelial cell specimen urine, sperm, blood serum/plasma	DNA	STI, HHV, HPV, Parvovirus B19, Borrelia and Anaplasma	yes
RealLine Extraction 100 CE-IVD	48	VBC8896	epithelial cell specimen liquor, biopsy material, blood serum/plasma, whole blood, urine, sperm, tick suspension	DNA and RNA	STI, HHV, HPV, Parvovirus B19, Toxoplasma gondii, Rubella, Borrelia, TBEV, Anaplasma HIV, HCV, HBV	yes
RealLine Extraction 100(A) For Research Use!	96	VBC8883	epithelial cell specimen liquor, biopsy material, blood serum/plasma, urine	DNA and RNA	STI, HHV, HPV, Parvovirus B19, Toxoplasma gondii, Rubella, Borrelia, TBEV, Anaplasma HIV, HCV, HBV	yes
RealLine Extraction 1000 For Research Use!	48	VBC8895	blood serum/plasma whole blood	DNA and RNA	HIV, HCV, HBV, Borrelia , TBEV	yes
RealLine FFPE DNA Extraction Kit For Research Use!	50	MED20301	FFPE tissue	DNA	Cancer Mutation Kits	no
RealLine FFPE DNA Extraction Mag For Research Use!	12	MED32501	FFPE tissue	DNA	Cancer Mutation Kits	yes

BIORON Diagnostics GmbH – Rheinhorststr. 18 - 67071 Ludwigshafen (Germany)

Phone: +49 621 545 900 70 Fax: +49 621 545 900 68

E-Mail: info@bioron.de Homepage: www.bioron.de